

EchoHeart® - the first transvaginal fetal Doppler probe in the world.

U.S. Patent #6379305

SummitDoppler
Listening to Life™

Peace of Mind in the Palm of Your Hand

EchoHeart helps streamline your practice by reducing the need for costly imaging exams when conventional Dopplers cannot detect the fetal heartbeat.

Prior to the EchoHeart, the only option was imaging.

Immediate benefits to your practice:


- Determine fetal heartbeats without scheduling an additional imaging exam
- Early detection, especially for obese or retroverted uterus patients
- Perfect for high-risk pregnancies when multiple exams are required

While the EchoHeart is not intended to replace imaging equipment, it can economically provide earlier fetal-heart detection for obese or retroverted uterus patients. In a recent study, clinicians using the EchoHeart were able to detect fetal heartbeats in 71% of patients between 6 and 8 weeks.

Since the probe shape and placement are similar to the transvaginal imaging probes used in most women's health care clinics, no special training is required for clinicians.

Special sheaths with preapplied coupling gel are available for patient safety.

The EchoHeart probe can quickly be interchanged with standard Summit Doppler 2 and 3 MHz probes, maximizing your equipment investment. (The EchoHeart will not work with Dopplers manufactured by other companies.)



For more information about our products, contact Summit Doppler Systems, Inc. at
4680 Table Mountain Dr. # 150 Golden, CO 80403 # 800-554-5090 www.summitdoppler.com

Customer service unequaled by any other Doppler manufacturer.

Day-to-day use can be hard on any equipment

Knowing that accidents with equipment can happen, Summit Doppler has created a customer service system to take care of the professionals that use our products. When you call, you talk to a live person – no endless phone menus to endure. If your Doppler ever needs to be repaired, we guarantee a 1-3 day turn-around on all repairs.

Full warranty

Additionally, all of our products are covered by a full, one-year warranty eliminating the "what-if" concerns that can accompany new equipment purchases. Your confidence and success is something we want to protect.

Common sense convenience for your office

Unlike many of our competitors, our products provide the convenience of interchangeable probes and user-replaceable cables. Offering a full line of probes, our Dopplers can provide increased efficiencies within any practice.

Summit Doppler

Summit Doppler Systems, Inc. is a team dedicated to Doppler ultrasound equipment. Doppler ultrasound is our flagship – not a back-burner product line. Together we have decades of medical equipment design and manufacturing experience. We are committed to using our extensive knowledge and experience to provide the finest products and service available. At Summit Doppler, we are determined to make the LifeDop your favorite Doppler – and Summit your Favorite Doppler Company.



The innovative EchoHeart transvaginal Doppler probe is the first of its kind – anywhere.