

Sterilizable Probe



SENSITIVE



RELIABLE



INTERCHANGEABLE

Summit Doppler Systems presents...

The LifeDop™ Sterilizable Probe for the detection of blood flow in both veins and arteries. This clinically tested probe is suitable for use in sterile fields. STERRAD® compatible for use with LifeDop Dopplers (model # L150, L150A, L150R, L250, L250R, L250AR & L350R only).

See other side for more information about the probe.

From the Doppler Manufacturer You Know and Trust...

The LifeDop[™] Sterilizable Probe

SummitDoppler
Listening to Life[™]



LifeDop 350 Table-Top Doppler with the sterilizable probe (model # STN)

**Contact Your Local Rep
or Summit Doppler for
more information.**

#800-554-5090

Summit Doppler
Products are
MADE IN THE U.S.A.

Key Features

- This sterilizable probe (model # STN) is an 8 MHz non-directional vascular Doppler transducer.
- This clinically tested probe with attached 7' straight cable is suitable for use in the sterile field.
- Sterilize using STERRAD[®] (an H₂O₂ gas plasma process) for up to 10 cycles.
- The STN probe interfaces with the Summit Doppler Systems LifeDop family of Dopplers, models: L150, L150A, L150R, L250, L250R, L250AR and L350R only.
- The sterilizable probe can be purchased individually or in packs of 3 for checking pulses and verifying distal blood flow.

Summit Doppler Systems, Inc.
4680 Table Mountain Drive # 150
Golden, Colorado 80403
800-554-5090 or # 303-423-7572
Fax # 303-940-7165

www.summitdoppler.com