

LifeDop™ 250 ABI (Model # L250AC)



Summit Doppler Systems offers fast, economical solutions for obtaining the ABI for diagnosing peripheral arterial disease (P.A.D.) that meet reimbursement criteria. The **LifeDop 250 ABI with 4 Cuffs & Aneroid** features a portable bi-directional Doppler with compact printer for the diagnosis of P.A.D. The system includes a hand-held Doppler, printer (waveforms required for reimbursement), 8 MHz bi-directional probe, adhesive-backed label paper, ABI chart, forms, ABI video, examination guide, quick reference guide and ABI glossary, carrying case, 4 blood pressure cuffs and aneroid package.

Features:

- Fast, waveform printout on adhesive-backed label paper for easy placement on reimbursement report forms
- Affordable, ultra-portable Doppler and printer with cable and recharger
- Display shows blood flow direction and relative frequency of the signal
- Automatic scaling and calibration
- Sensitive, bi-directional probe and hand-held Doppler
- Dimensions (h w l): 140 x 70 x 35 mm; Weight: 320 grams
- Full 1 year product warranty

Summit Doppler Systems, Inc.
4680 Table Mountain Dr. #150
Golden, Colorado 80403
303-423-7572 or 800-554-5090
Fax # 303-940-7165
www.summitdoppler.com